

Bruker AVIII HD 400-C3PO and 500-Chewie Training Note:

1. Open topspin
2. Start → Create Dataset (new)
 - Select solvent; check 'getprosol'; select pulse sequence: PROTON or C13CPD
3. Acquire →
 - Sample (ej & ij) --- clean nmr tube and hold spinner properly.
 - Lock (lock) --- wait till "lock done"
 - Tune (atma)
 - Spin (ro on; ro off; ro)
 - Shim (topshim)
 - Prosol (getprosol)
 - Gain (rga)
 - Then you may go to "AcquPars tab" change 'NS'; and click "clock" icon to show total experiment time
 - Go (zg)
4. Process →
 - Proc. Spectrum ('procl d y' or efp)
 - Adjust Phase (apk)
 - Calib. Axis (.cal)
 - Pick peaks (.pp)
 - Integrate (.int)
 - Correct baseline (absn)

Note:

- 1) Always **PATIENTLY** wait till "tune or other process done or finished" shown up and then go to next step! Important!
- 2) Drag previous data folder in, and "use current parameters"...; always go through the above procedure once again for every new sample; must do **TUNE** if you change to different nucleus 1D measurement, e.g. from PROTON to C13CPD.
- 3) Upload your data from (/home/your_nmr_username/nmr --- on B400-C3PO; or /opt/topspin3.5plx/data/your_nmr_username/nmr --- on all other Bruker machines) into your OneDrive.
- 4) If you want to run other 2D homonuclear (COSY, NOESY, ROESY, TOCSY, etc.) or heteronuclear experiments (HSQC, HMBC, etc.), ask Hongwei first before you run.
- 5) A more detailed instruction of Bruker topspin can be found on our "GT nmr website" under "forms and instruction" tab.