

## Bio-NMR user Training Procedure on Bruker AVIII HD 800 and 700:

1. Open topspin
2. Start → Create Dataset (new)
  - Check Experiment and select ZGPR; Select solvent (H<sub>2</sub>O+D<sub>2</sub>O) for 90%H<sub>2</sub>O/10%D<sub>2</sub>O buffer; check 'getprosol'.
3. Acquire →
  - Sample (ej & ij) --- wipe and clean nmr tube and handle spinner and sample gage properly. (sx 3 or sx 16 for the sample in # 3 or # 16 position.)
  - Lock (lock) --- wait till "lockn: done."
  - Tune manually (atmm) for samples using 90%H<sub>2</sub>O/10%D<sub>2</sub>O. (atma) --- automatically tuning.
  - Spin (ro off; ro on; ro); don't spin for protein samples.
  - Shim (topshim; topshim tune tunea)
  - Prosol (getprosol);
  - for ZGPR, use (pulsecal) to calibrate 1H 90° pulse -- P1; set rg = 1 and (gs) to find O1, which is the center of the spectrum, and water peak position too.
  - Gain (rga) — automatically rg adjustment.
  - Go to "AcquPars tab" to change 'NS'; and click "clock" to show total experiment time.
  - Go (zg);
  - load ZGESGP, update P1, O1, and set proper NS, and then zg.
  - HSQCETFPF3GPSI is a good one for 2D N15 HSQC...
4. Process →
  - Proc. Spectrum ('proc1d y' or efp)
  - Adjust Phase (apk)
  - for 2D: use (xfb – 2D FT); (rser 1 – read first FID) (efp, apk, .ph)
  - Calib. Axis (.cal)
  - Pick peaks (.pp); Integrate (.int); Correct baseline (absn)

### Note:

- 1) Always **PATIENTLY** wait till "tune or other process done or finished" shown up and then go to next step! Important!
- 2) must do **TUNE** again if you use other nuclei (C13 /N15) besides H1.
- 3) always go through the above procedure once again for every new sample.
- 4) Drag previous data folder in, and "use current parameters" ...
- 5) Upload your data from (/opt/topspin3.5pl5/data/your\_nmr\_username/nmr --- on B800 and B700) into your OneDrive.
- 6) Always ask Hongwei first before you run any new things.